

2012 Influenza Results Interpretation Michigan Department of Community Health, Bureau of Laboratories

Report Appears As:	Interpretation
<hr/> TEST RESULTS <hr/> <p>COMPREHENSIVE INFLUENZA PCR PANEL INFLUENZA 2009 A/H1 RNA DETECTED</p> <p>Testing performed by Real Time RT-PCR.</p> <p>Fact sheets for health care providers and patients are available at www.cdc.gov/h1n1flu/eua/testkit.htm.</p>	<p>The patient has the Influenza 2009 A/H1 virus. This virus was previously known as the 2009 Novel Influenza A (H1N1) virus. For Influenza 2009 A/H1 RNA to be detected, all influenza tested markers; Universal Influenza A, Pandemic Influenza A, Pandemic H1 and internal control are amplified.</p>
<hr/> TEST RESULTS <hr/> <p>COMPREHENSIVE INFLUENZA PCR PANEL INFLUENZA A AND B RNA NOT DETECTED</p> <p>A negative result does not preclude infection by another respiratory pathogen.</p> <p>Testing performed by Real Time RT-PCR.</p> <p>Viral culture results to follow.</p>	<p>No indication of any Influenza virus. The specimen has been placed in viral culture to test for other respiratory viruses.</p>

<p style="text-align: center;">TEST RESULTS</p> <hr/> <p>COMPREHENSIVE INFLUENZA PCR PANEL INFLUENZA B RNA DETECTED</p> <p style="text-align: center;">Testing performed by Real Time RT-PCR.</p>	<p>The patient has Influenza B virus. This is a seasonal Influenza virus.</p>
<p style="text-align: center;">TEST RESULTS</p> <hr/> <p>COMPREHENSIVE INFLUENZA PCR PANEL SEASONAL INFLUENZA A, SUBTYPE H3 RNA DETECTED</p> <p style="text-align: center;">Testing performed by Real Time RT-PCR.</p>	<p>The patient has Influenza A subtype H3 virus. This is a seasonal Influenza virus.</p>
<p style="text-align: center;">TEST RESULTS</p> <hr/> <p>COMPREHENSIVE INFLUENZA PCR PANEL SEASONAL INFLUENZA A, SUBTYPE H1 RNA DETECTED</p> <p style="text-align: center;">Testing performed by Real Time RT-PCR.</p>	<p>The patient has Influenza A subtype H1 virus. This is a seasonal Influenza virus.</p>

<p style="text-align: center;">TEST RESULTS</p> <hr/> <p>COMPREHENSIVE INFLUENZA PCR PANEL TEST NOT DONE</p> <p>Quantity insufficient for testing. This specimen source has not been validated for this assay.</p>	<p>Specimen received was either leaking or did not meet testing guidelines.</p>
<p style="text-align: center;">TEST RESULTS</p> <hr/> <p>COMPREHENSIVE INFLUENZA PCR PANEL INCONCLUSIVE</p> <p>Specimen sent to CDC for further testing. Subtype could not be determined by PCR.</p> <p>Testing performed by Real Time RT-PCR.</p>	<p>Comprehensive Influenza virus PCR indicated presence of Influenza A, but failed to subtype.</p>
<p style="text-align: center;">TEST RESULTS</p> <hr/> <p>COMPREHENSIVE INFLUENZA PCR PANEL INFLUENZA A RNA DETECTED, PRESUMPTIVE POSITIVE FOR H3 VARIANT</p> <p>Specimen sent to CDC for confirmation. Testing performed by Real Time RT-PCR.</p>	<p>Comprehensive Influenza virus PCR indicated presence of universal and pandemic Influenza A, presence of H3 but failed to indicate the presence of seasonal and pandemic H1 subtype</p>

TEST RESULTS	
<p>COMPREHENSIVE INFLUENZA PCR PANEL</p> <p style="text-align: center;">EQUIVOCAL</p> <p style="text-align: center;">Specimen sent to CDC for further testing.</p> <p style="text-align: center;">Subtype could not be determined by PCR due to low viral titer.</p> <p>Testing performed by Real Time RT-PCR.</p>	<p>Influenza virus markers failed to amplify according to test protocol such that a result could not be determined.</p>

Results with double asterisks indicate a test was reflexed off this result OR the result was unexpected in this assay, e.g., equivocal or inconclusive.

For all inquiries, contact Anthony Muyombwe, PhD Virology Section Manager MDCH, Bureau of Laboratories 3350 N MLK Jr. Blvd Lansing, MI 48909 Tel: 517-335-8099 Fax: 517-335-9631 Email: muyombwea@michigan.gov